


freezeTb; detecting anti-microbial resistance in TB


WHO 2023



+





=



tuberculosis

5-7 million 10 million 1 million


Treatment success



WHO 2023


Drug resistant detection


Month





Campelo et al. 2021

Campelo et al. 2021;
Mukherjee et al. 2023








Does not include van or great danes


Xia et al. 2023

freezeTb; the software side

Expected output



Check for AMRs




Detect Lineages

2. a. MIRU-VNTR

- Amplicons target 24 repeat regions across the genome
- Amplicons lineage found by Length of amplicon
- Current gold standard

2. b. Spoligotyping

- Amplicon targets direct repeat region (1 primer)
- Lineage found by detecting 43 internal spaces
- Lineage id is barcode with presence/absence of spacer



Build consensus

TBProfiler; the current resource

TBProfiler; pros

- Detects AMRs uses WHO's 2023 catalog
- Detects spoligotypes and SNP lineages (Nanopore?)
- Outputs an vcf files

TBProfiler; cons

- More focused on Illumina
- > 5 dependencies or conda; not fun to install

Other solutions;

ONTs: Have not tried

- Detects spoligotypes
- Likely more to add

Plenty of small programs that add in more dependencies

AMR detection; no small program I know of

Lineage detection is simple


Consensus building; non-model methods are simple

Most of it is simple, why not make our own?


FreezeTb; the software side

- CLI; requires sam file
- R GUI; fastq file only, but requires minimap2

Pre analysis




Mapped reads




Length < 200
Q-score < 7
Mapq < 15


Filter reads



Read depth



Coverage



Build consensus

Filters out bases with Q-scores under 7

Break up consensus

- Fragments assemblies under 20x read depth
- Remove fragments less than 200 bases long


Masking:

- Select nucleotides with most support
- Mask nucleotides with less than 50% of mapped reads

Indel removal:

- Remove indels with less than 70% of mapped reads

Output and sam file



Check for AMRs

Using 2023 WHO tuberculosis AMR mutation catalog

Gene deletions are ignored

Frames shifts / LoF's treated as normal AMR


- Less false positives due to indel errors

For amino acid catalog variants

- AMR region in read converted to amino acids

AMR called if

- Sequences match
- Read AMR region length = catalog AMR length



Detect Lineages


Miru lineages:

- See if read has primer coordinates (from sam file)
- Check reads repeat length

Spoligotyping:

- Get direct repeat (DR) region from consensus
- Map spacer sequences to DR with Waterman alignment
- Make barcode from mapped spacers and find lineage

Pre analysis




Mapped reads


Length < 200
Q-score < 7
Mapq < 15
Filtler reads



Build consensus



Read depth Coverage



Detect Lineages

Check for AMRs

Each step in freezeTb is modular:

- You like an module; you can get an module
- Each module takes an sam file

FreezeTb; works, but not an replacement

Future steps:

- Finish and integrate spoligotyping
- Maybe add in fast kmer seeding for spoligotyping

Fast kmer seeding is close to an read mapper

Why not add an read mapper? Ok probably not.

